

## Drive DeVilbiss 10L Oxygen Concentrator

HCPCS: E1390



FIGURE A



FIGURE B



FIGURE C

- Built upon the reliable features of the 525 concentrator, the 1025 Oxygen Concentrator is designed to provide optimal oxygen delivery for patients with varying flow requirements
- Assembled in the USA and built in the same shell as the 5 Liter, the 10 Liter Concentrator is also compatible with cylinder transfill systems
- Sold separately, the Transfill Caddy is designed with a steel frame to simply connect the 18 x 14 inch shelf, designed to support up to 40 lbs., to the DeVilbiss 10L Concentrator
- Front view (Figure A), Rear View (Figure B)
- With a weight of 42 lbs. and an identical footprint to the 525 unit, the 1025 is substantially smaller and lighter than competitive 10 Liter concentrators
- Capable of delivering 87% to 96% of O<sub>2</sub> purity from 2 to 10 LPM, the 10 Liter Concentrator is an ideal investment for your patients and business alike
- Exclusive DeVilbiss OSD<sup>®</sup> (Oxygen Sensing Device), standard on every unit, monitors the oxygen produced and ensures patient safety and reliability for longer service intervals
- Equipped with accessible patient controls, bright LEDs that can be seen from a distance, protected cannula fitting and recessed humidifier nook to prevent damage

### PRODUCTS

#### Description

1025DS	Drive DeVilbiss 10L Oxygen Concentrator, EA
--------	---

### SPECIFICATIONS

#### Description

Delivery Rate	2.0 to 10.0 LPM
Outlet Pressure	20.0 PSI $\pm$ 1.0 PSI
Oxygen Percentage	2-10 LPM, 87-96%
Power Consumption	120VAC, 60Hz
Power Input	639W (avg)
Product Length	24.5"
Product Overall Height (FW17   FW19)	12"
Product Weight 2   4   6	42 lbs.
Pulse Rate	13.5"
Warranty Limited	3 Years

### ACCESSORIES

#### Description

525DD-650	Transfill Rack System, 525DS Series, 1/e, EA
PF1100TUB	Transfill Connection Tubing, 1/ea, EA
MQ0412	No Smoking Sign, 25/CV, CS
HUM 001	Humidifier Bottle PI 50/CS, CS

**PARTS****Description**

303DZ-605	AIR FILTER 6-PKG - DFT, PKG
PV5LD-651	FINAL BACTERIA FILTER - DFT, EA
1025D-605	INTAKE FILTER SSP, 1/ea, EA
1025D-682	COMPRESSOR FILTER SSP, 1/ea, EA