



Fall Management

Supply Staff with Effective and Cost-Efficient Fall Detection

Sensatec solutions keep people safe, improve care, and increase staff efficiency. Sensatec pads and monitors are a simple fall prevention solution that helps staff provide the highest level of care.



Easy-to-use Automated Operation


Start using the product right away. SENSATEC® is easy to operate and involves minimal staff training.

The auto-on feature automatically begins monitoring immediately after a patient is positioned on a bed or chair sensor pad.

The auto-reset feature intuitively resets the unit when a patient's weight is reapplied to the pad.



Advanced Technology Automatically Alerts Potential Falls



Comfortable, non-slip sensor pads are available in multiple models for short in-patient stays or longer-term monitoring.

ST750 Fall Monitor Features & Benefits

- Audible alarm options: choose from eight melodies, eight tones, or a 10-second voice recording
- Alarm delay can be selected for zero, one, two, or four seconds
- Adjustable volume for different proximity to staff: Low, Medium, High, or Silent
- Nurse call integration
- Optional AC adapter
- Alarm units undergo rigorous quality testing and have a 1-year warranty

