

SENSATEC®

Fall Monitors

Begins monitoring immediately when a patient is positioned on a sensor pad, triggers an alert when weight is removed, and intuitively resets when a patient's weight is reapplied. Can easily hang or hook on beds, chairs, or wheelchairs.



Dimensions (LxWxD)

4.7" x 2.8" x 1.3"

Weight

0.4 lbs

Power Requirement

(3) AA Alkaline Batteries
Power Adapter, 5VDC, 500mA (optional - ST750 only)

Jack Connections

RJ9 Modular Sensor Input Jack
0.25" Stereo Floor Sensor Jack
0.25" Stereo Nurse Call Jack (ST750 only)

Audio Output

78dBm

Temperature

Operating: 41-104°F, 95% non-condensing

Warranty

One year

Status Indicator

Solid Green: 30 second Hold
Flashing Red: In Alarm
Blinking Red (every 5 seconds): Low Battery

Sensor Pads

Oversized sensor pads ensure reliable monitoring and put an end to false notifications. Durable latex-free pads are constructed of an innovative non-slip material that adheres to a mattress or chair without the need to peel-off adhesive strips or rubber bands.

Dimensions (LxW)

30" x 15"

Weight

7.3 - 10.2 oz (depending on life of pad)

Composition

Polyvinyl

Cord Length

54" (connects to fall monitor)

Dimensions (LxW)

15" x 12"

Weight



3.5 - 4.2 oz (depending on life of pad)

Composition

Polyvinyl

Cord Length

24" (connects to fall monitor)



Dimensions (LxW)

30" x 15" (bed) or 15" x 9" (chair)

Weight

0.7 - 2.2 lbs

Composition

Antibacterial Polyvinyl

Cord Length

Bed: 72" (connects to fall monitor)
Chair: 42" (connects to fall monitor)



Dimensions (LxW)

36.6" x 24.6"

Weight

3.95 lbs

Composition

Polyvinyl

Cord Length

103" (connects to fall monitor)

SENSATEC®

Additional Sensor Options

Seatbelt sensor and toilet sensor pad options are available for greater protection.

Dimensions (LxW)	52" x 2"
Weight	6.2 oz
Composition	Latex free, nylon, velcro
Cord Length	40" (connects to fall monitor)
Storage Temperature	50-104°F, 30-75% relative humidity



Dimensions (LxW)	5.3" x 2.8"
Weight	1.3 oz
Composition	Polyvinyl
Cord Length	73" (connects to fall monitor)



Accessories



Mounting Clip



Wall Bracket



Nylon Strap



Velcro Strap & Strip



10' Extension Cord



Nurse Call Y Adapter



Replacement Boot



AC Adapter



Nurse Call Connection Cord

