



Crest GeriCall™ Call Cords



Description

The GeriCall Geriatric Call Cord is a round, flattened momentary call cord pad that works great for residents who have limited hand dexterity. This call cord utilizes a mechanical switch in the pad that can be activated by slight pressure from the hand, arm or body, or can even be positioned under the resident's chin.

Features

- Large pad for easy activation
- Includes security clip
- 3" diameter pad
- FDA Listed

Technical Specifications

Electrical Rating:	28 VAC/DC	Color:	Classic Grey
Electrical Connection:	Various connection types available	Weight:	Varies based on cord length
Switch Type:	Normally Open SPST	Cord Length:	8' - 15'
Warranty:	2 Years	Dimensions:	1 ⁵ / ₃₂ "(H) x 3 ³ / ₁₆ "(W) x 4 ¹³ / ₆₄ "(L)
Material:	Silicone Synthetic Rubber	Environment:	Operating temperature: 50° - 100°F
Standards:	RoHS Compliant FDA Listed		Relative humidity: 30% to 80% (non-condensing)

Crest Technical Support: **1-800-328-8908**

Specifications subject to change without notice.
Visit CrestHealthcare.com for additional technical documents.