



Crest DuraGlo™ Momentary Call Cords



Description

DuraGlo is the only call cord in the industry with glowing bands that help your residents locate their call cord in dark rooms! DuraGlo requires no power to glow - it works with any nurse call system.

DuraGlo can glow for up to 4 hours in total darkness, depending upon the light source used and the length of time exposed, and will recharge with every exposure to a light source. Best of all, DuraGlo is built on the same solid platform as our popular DuraCall call cords.

Features

- Bright glow requires no power
- Only call cord with glowing bands
- Locate call cord easier in dark rooms

Technical Specifications

| | | | |
|-------------------------------|--------------------------------------|---------------------------|---|
| Electrical Rating: | 24 VAC/DC | Weight: | 5.8 oz (7' version) |
| Electrical Connection: | ¼" Plug | Pendant Dimension: | 2.75"H x 1.625"W |
| Button Life: | >1,000,000 cycles | Cord Length: | Various lengths available |
| Switch Type: | Normally Open SPST | Environment: | Operating temperature: 50° - 100°F |
| Warranty: | 2 Years | | Relative humidity: 30% to 80% (non-condensing) |
| Material: | UL 94V-0 Flame Rated, Latex Free PVC | | Oxygen safe: Meets UL 1069 oxygen enriched atmosphere testing |
| Standards: | FDA Listed RoHS Compliant | | Liquid ingress protection: Meets IEC 60529 IPx4 requirements |
| Color: | White pendant with red call button | | |

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit CrestHealthcare.com for additional technical documents.