

Expander

Interfaces fall management into existing nurse call system



The Expander connects wired and wireless sensing devices to your existing nurse call system(NCS).

MEX-10 is a versatile but economic NCS splitter with disconnect alert. It allows user to connect up to 3 wired accessories into a single call point. It is designed to accommodate the most common NCS in hospital & elderly care facilities, includes floor mat input.

CMEX-11 is a cordless version. It allows the use of up to 5 Rondish wireless devices in the patient room. This product also allows a hospital to bypass a standard bed monitor and output exit alarms directly to their NCS. Eliminating local alarms helps to promote a peaceful recovery environment.

CMEX-21 provides the same benefits with CMEX-11, while allowing flexibility in connection types. In some cases 6.3mm (1/4") plug connection may not fit the existing NCS socket. Therefore the CMEX-21 comes with a separate cable with plug, allows the flexibility to change the connector types used by nurse call providers.

Extra features of Expander

MEX-10 (Expander)



- 2 x 6.3mm (1/4") mono mic jack ports for call cord
- Use call cord cable to connect with bed monitor
- 1 x RJ10 port for floor mat
- Allows user to plug 3 accessories into a single call point
- Cord removal alert (with on/off switch to disable this function for Port 2)
- Compatible with normal open system
- Compatible with mono/stereo type nurse call panels
- No power required

CMEX-11 (Cordless Expander)



- 433 MHz receiver allows use of floor mat, bed / chair pad, nurse call handset, wristband transmitter, bathroom alarm, motion sensor & door alarm, up to 5 wireless devices
- 1 x 6.3mm (1/4") mono mic jack port for call cord
- 1 x RJ10 port for wired floor mat or bed pad
- Cord removal alert (with on/off switch to disable this function)
- Compatible with normal open or normal closed systems
- LEDs to indicate the alarm status & error types, such as lost signal, out of range and low battery indication
- Button on the top for clearing alarm or pause function under no alarm condition
- Pause function allows patient removal from bed
- 3 x AAA batteries or 5V DC adapter

CMEX-21 (Cordless Expander with Interchangeable Plug)



- Interchangeable cable plug allows the flexibility for different connector types nurse call systems (Cable can be made to order)
- 433 MHz receiver allows use of floor mat, bed / chair pad, wristband transmitter, bathroom alarm & door alarm, up to 5 wireless devices
- 1 x 6.3mm (1/4") mono mic jack port for call cord or bed monitor
- Cord removal alert (with on/off switch to disable this function)
- Compatible with normal open or normal closed systems
- LEDs to indicate the alarm status & error types, such as lost signal, out of range and low battery indication
- Button on the top for clearing alarm or pause function under no alarm condition
- Pause function allows patient removal from bed
- 3 x AAA batteries or 5V DC adapter
- A hanger clip & a wall mount bracket is included

