

***Product*** Receiver Dongle

**Product Code:** DON-30

**Description:** Receiver dongle for Nexus. Connects Nexus with LoRa wireless devices for rapid deployment of a basic nurse call system



**Features:**

- Outputs signal data to USB-type A port
- Includes software license for Nexus
- Handshakes with indication light to confirm signal is received

**Works with:**

- NDP-11 Nexus Display Panel
- TVB-21 Nexus Control Box
- ILB-21 Indication Light

***Technical Specifications***

Operating frequency: 435MHz  
Operating range: Up to 1,500m (line of sight)  
Power supply: Powered from tablet  
  
Dimensions: 8 x 2 x 1cm  
Housing material: ABS  
Net weight: 10g  
Operating temp.: 5°C to 40°C  
Rel. humidity: <90%

***Standards***

EMC: EN 55022  
Radio: EN 300220-1&2,  
EN 301489-1&3  
Alarm: EN 50130-4  
Safety: EN 60950-1  
  
Manufacture: ISO 9001:2008,  
ISO 13485:2003