



Pullcord Station, Crest Replacement for Couch™

HA-1334-1C

Typical Installation

- White — (L2) Lamp return common, system plus
- Blue — (D) Direct wire to annunciator
- Red — (CL) Direct wire to corridor light
- Yellow — (L1) Normal call system common

Specifications

ELECTRICAL

Current & Voltage Rating: Resistive load 1.5Ap @ 28VDC.
0.50A @ 115VAC.

MECHANICAL

Terminals: Brass, silver plated
Contact Life: 25,000 cycles minimum
Mounting: Fits standard single gang back box

SIZE MEASUREMENTS

Plate: 4.5" x 2.75" (1 Gang)

Function Table

Color	Violet	Grey	Red	Yellow	White	Blue
Original Terminal Designation	3	5	4	6	1	2
Up Position	On	On	Off		Off	Off
Middle Position	Off	Off	On		On	On
Down Position (Momentary)	Off	Off	On		On	On
Green wire is ground						

Electrical Schematic

