

Features

- Easy Installation
- Mounted on a flame-retardant (UL® 94V-0) plastic panel
- Simple, reliable operation

Description

The JK300A Monitor Call Interface is for use in applications where a simple interface is needed to connect the monitor input of an IR300-series or IR400P5-series station to an external contact closure. The JK300A accepts a standard ¼" phono-type jacks and places a call if a plug is removed. A ¼" dummy plug is included for use when the monitor point is not connected to an external device. The JK300/2 Monitor Call Interface provides two ¼" jacks and two dummy plugs.

Specifications

Dimensions:	Height:	4.5"	(114 mm)
	Width:	2.75"	(70 mm)
	Depth:	1.5"	(38 mm)
	Projects:	0.25"	(6.5 mm) from mounting surface
Construction:	Reliable components mounted on a flame-retardant (UL® 94V-0) plastic panel		
Connections:	Wire lead pigtailed		
Housing:	Single-gang electrical box or ring		



JK300A, JK300/2 Monitor Call Interface

Applicable Systems

Tek-CARE®400 P5/P5+ Nurse Call System
IR400P5/P5+-series Multipurpose, Patient, Resident, Check-In, Staff or Duty Stations
-or-
Tek-CARE®400 Nurse Call System
PM423 Auxiliary Input Module
IR400-series Multipurpose, Patient, Resident, Check-In, Staff or Duty Stations
-or-
Tek-CARE®NC300III Nurse Call System
IR300-series Multipurpose, Patient, Staff or Duty Stations
-or-
Tek-CARE®120 or 160 Nurse Call Systems
PM123 Auxiliary Input Module