



VIVE
TRANSPORT WHEELCHAIR
MOB1021BLK
HCPCS: E1038

The Vive Foldable Transport Wheelchair safely supports up to 300 pounds, with a wide, 19" seat and padded back and armrests for exceptional comfort. The transport chair features two folding points to fold flat for compact storage and easy transportation in most vehicles. Constructed with durable aluminum, the wheelchair includes 12-inch tires and dual hand brakes for greater stability and safety on any surface.

WHAT'S INCLUDED

- Transport Wheelchair
- Lifetime Guarantee

COLORS   

FOLDABLE TRANSPORT WHEELCHAIR



Specifications

Base when open	21" by 32"
Handle Height	40
Front Wheels	8"
Back Wheels	12"
Wheels Diameter	8"
Seat Height	19.5"
Folded	8.25" by 30"
Seat Widths	19"
Seat Depth	16"
Weight limit	300 pounds



SAFE, COMFORTABLE TRANSPORTATION

Providing safe, comfortable mobility for adult patients and loved ones, the Vive foldable wheelchair features a wide 19" seat, padded armrests, a quick release seat belt and a breathable, padded backrest. The transport chair is easy to maneuver both indoors and out with large, 12-inch rear tires for greater stability and traction.

EASILY FOLDS FLAT

The transport wheelchair folds flat for convenient storage and travel. Compact for easy portability, the backrest and handle section also fold down to easily fit in most vehicles.

DURABLE ADJUSTABLE FRAME

Pre-assembled, the wheelchair frame is constructed with durable aluminum that safely supports up to 300 pounds. The foldable wheelchair features adjustable footrests that can easily be raised or lowered for greater patient comfort.

ERGONOMIC HANDLES

Designed to reduce hand and wrist fatigue, the handles are gently angled outward and are covered with a textured grip for greater comfort. The transport chair also includes loop hand brakes to lock the wheels during transfers.

VIVE GUARANTEE

Lifetime guarantee so you can purchase now with confidence.

