



**PROVIDING COMFORT
AND STABILITY
WHILE WALKING**

VIVE
UPRIGHT WALKER
MOB1033

The Vive upright walker features tall adjustable handles to eliminate hunched shoulders and promote proper posture. Each arm rail is height and length adjustable and features soft, contoured armrests and textured hand grips for greater control and comfort. The lightweight aluminum walker weighs approximately 20 pounds and is compact when folded for easy travel.

WHAT'S INCLUDED

- Upright Rollator Walker
- Removable shoulder bag
- Cane Holder
- One year guarantee

COLORS





portal.vivehealth.com



SUPPORTS WALKING UPRIGHT

Eliminating tension in the shoulders and neck, the Vive upright walker is height adjustable, allowing individuals to stand up straight and look forward while walking to promote a natural gait. The upright walker rollator features 8", flat-free wheels for greater stability on moderately smooth surfaces.

EASY HEIGHT ADJUSTMENTS

Locking push pins allow the height of the Vive upright walker handles to be easily adjusted from 45" to 52" in height. The folding walker handles are topped with contoured armrests to evenly distribute pressure, relieving fatigue and length adjustable hand grips with easy to use loop brakes to lock the wheels in place for stable seating.

COMFORTABLE ON THE GO SEATING

Featuring a supportive, 17.5" wide seat, the tall upright walker provides comfortable, on-the-go seating for resting as needed. The upright walker also includes a detachable shoulder bag that safely stores essential items close at hand.

LIGHTWEIGHT AND FOLDABLE

Constructed with a lightweight aluminum frame, the Vive upright walker is exceptionally durable, safely supporting up to 300 pounds. Easy to transport, it folds for compact storage, fitting in the trunk of a car or even in the backseat.

Specifications

Overall Width	25"
Overall Length	32"
Armrest Lowest Height	41"
Armrest Highest Height	46"
Seat Width	17.5"
Seat Depth	9.5"
Folded Width	10.8"
Removable Bag Weight Capacity	11 lbs.
Product Weight	20 lbs.
Max User Weight	300 lbs.

vive health 8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

SR22- 6R1



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com