

Product Nexus Rugged Tablet

Product Code: NDP-21

Description: Rugged tablet with Nexus software for basic call system. Ideal for rapid deploy nurse call. Must include receiver dongle or ZP-01(M) gateway to operate



Features:

- 10" Industrial-grade Android tablet with docking station for convenient placement on nurse desk
- Receiving options for smaller, medium and large wireless systems
- Smartphone app available with acknowledgement
- Signaling options for remote displays or pagers (requires docking station)

Works with:

- DON-20/30/40 Wireless Receiver Dongle
- ZP-01(M) Master Protektor Unit
- GP2009TR Pager Controller
- DOCK-01 Nexus Docking Station

Technical Specifications

Radio frequency: Varies according to dongle
Operating range: Varies according to dongle
Power supply: DC5V input, 2000mA
Battery back-up for 8+ hours
Dimensions: 280 x 190 x 25mm
Net weight: 1250g
Operating temp.: 5°C to 40°C
Rel. humidity: <90%

Standards

EMC: EN 301 489-1, 3 &17
EN 55022,
EN 55024
EN 50332-2
EN 62479
EN 300 328
Radio: EN 300 440-1&2
Safety: EN 60950-1