

**Product** Nurse Call Handset

**Product Code:** NEC-12M

**Description:** Nurse call handset with call button and wall-mount bracket. Connects with Bedhead Call Point

**Features:**

- Wipeable surface for easy sterilization
- Reassurance LED
- Plastic wall-mount bracket
- USB break point to protect products

**Works with:**

- WCP-11/21C Bedhead Call Point



**Technical Specifications**

Dimensions: 12 x 5 x 2.5cm (handset only)  
Housing material: ABS  
Net weight: 160g (with cable and bracket)  
Operating temp.: 5°C to 40°C  
Rel. humidity: 85%