

Product Safety First Monitor

Product Code: NGM-41

Description: Bed exit monitor alerts caregiver to patient activity; wireless link to sensorpads, pendant, floor mat, or motion sensor

Features:

- Works with up to 3 wireless devices
- Pause, reset & volume control on monitor face
- Convenient On/Off switch on left side
- 6.3mm audio output to trigger nurse call
- NO/NC switch in battery compartment
- Operates on batteries or AC adaptor
- Low battery and change pad notification
- Input lost notification (pad lost, nurse call)

Works with:

- PT-11 Snap-on Sensorpad Transmitter
- TXP-11 Emergency Transmitter w/ Man-Down
- PIR-21 Wireless Motion Sensor
- MDC-11 Magnetic Door Sensor
- PC-40 Protective Cover
- BKT-10 Wall Mount Bracket



Technical Specifications

Radio frequency: 433.92 Mhz
 Operating range: 80m (line of sight)
 Power supply: 5~6V DC, 300mA
 or 3 x AA-type alkaline batteries
 Battery life: 6 months (Nominal)
 Dimensions: 120 x 83 x 35mm (w/ cover)
 Housing material: ABS (UL-94HB)
 Net weight: 210g (w/ cover)
 Operating temp.: 5°C ~ 40°C
 Rel. humidity: < 85%

Standards

EMC: EN 55022
 Radio: EN 300220-1&2,
 EN 301489-1&3
 Social Alarm: EN 50130-4
 Safety: EN 60950-1
 Medical Device
 Immunity: EN 60601-1-2
 Manufacture: ISO 9001:2015,
 ISO 13485:2016