

Product Emergency Transmitter w/ Man-Down

Product Code: TXP-11

Description: Worn as either wristband or pendant, and allow patients to call for emergency alarms. Programmable automatic fall-alert is available

Features:

- Active emergency button sends signal to central monitor or field receiver
- Includes lanyard with anti-ligature point
- Optional watch strap (not included)
- Accelerometer to trigger call in case of patient fall (recommended for pendant only)
- Reset signal triggered using magnetic switch
- Changeable battery - uses CR2450
- Waterproof to IP65
- Optional RFID tag for access control

Works with:

- ZP-01(S) Zone Protektor
- CMEX-11/21 Cordless Expander
- NGM-21 CordFree Monitor
- PAG-11/21/31 Caregiver Alert Pager
- ILB-21 Indication Light Box
- CK-01 Caregiver Key



Technical Specifications

Operating frequency: 433.92MHz
Transmission power: 10dBm
Operating range: 50-100m (line of sight)
Power supply: 3Vdc, Li button cell, CR2450
Battery life: 1-2 years
Dimensions: 46 x 34 x 13mm
Housing material: ABS + PC
Net weight: 25g
Operating temp.: 5°C to 40°C
Rel. humidity: 100%

Standards

EMC: EN 55032
Radio: EN 300220-1&2,
EN 301489-1&3
Alarm: EN 50130-4
Safety: EN 60950-1
Medical Device
Immunity: EN 60601-1-2
Manufacture: ISO 9001:2008,
ISO 13485:2003
IP rating 65
FCC Approved WNG-TXP-11